



The NEW Synea Fusion: A synthesis of design, technology and value

NEW!



The new range of Synea Fusion handpieces from W&H builds on the long established and respected reputation of the W&H Synea brand. Synea Fusion offers a practical range of quality handpieces at an exceptionally attractive price.

Love handpieces? Love Synea Fusion.

syneac FUSION



Type:	TG-98 L / TG-98	TG-97 L / TG-97	TG-98 LM	TG-97 LM
Light:	LED + / no	LED + / no	coupling dependent	coupling dependent
Coupling:	Roto Quick	Roto Quick	Multiflex*	Multiflex*
Head size (Ø):	11.5	10.0	11.5	10.0
Power (W):	19	16	19	16
Speed (rpm):	360,000	390,000	360,000	360,000
Max . length of rotating instruments (mm):	25	21	25	21
Operating pressure (bar):	3 ± 0.3	3 ± 0.3	2.5 – 4	2.5 – 4
Spray:	4-hole	4-hole	4-hole	4-hole

* Multiflex® is a registered trademark of Kaltenbach & Voigt GmbH, Germany



Type:	WG-99 LT / WG-99 A	WG-66 LT / WG-66 A	WG-56 LT / WG-56 A	HG-43 A
Light:	halogen / no	halogen / no	halogen / no	no
Head size (Ø):	10.2	9.5	9.5	–
Transmission/reduction:	1:5	2:1	1:1	1:1
Maximum motor speed (rpm):	40,000	40,000	40,000	40,000
Rotating instruments:	FG-burs	contra-angle burs	contra-angle burs	handpiece and contra-angle burs
Motor coupling:	ISO 38964	ISO 38964	ISO 38964	ISO 38964
Spray:	4-hole (> 50ml/min)	1-hole (> 50ml/min)	1-hole (> 50ml/min)	1-hole (> 50ml/min)
Intensity of illumination (light versions):	25,000 lux	25,000 lux	25,000 lux	–

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
 Ignaz-Glaser-Straße 53, Postfach 1
 5111 Bürmoos, Austria
 t + 43 6274 6236-0
 f + 43 6274 6236-55
 office@wh.com
 wh.com

Distribution:

W&H (UK) Limited
 6 Stroud Wood Business Centre
 Park Street, St Albans, Hertfordshire AL2 2NJ, United Kingdom
 t + 44 1727 874990
 f + 44 1727 874628
 office.uk@wh.com
 wh.com

