




**NEW**  
**DÜRR VISTASCAN**  
**MINI VIEW**  
**2.0**



CBCT Imaging  
Surgical  
X-Ray Equipment  
Dental Chairs  
Plant Equipment



**EQUIPMENT**  
**SALES**  
**& SERVICE**



X-Ray Compliance  
Scheduled  
Maintenance  
Equipment Repairs




\*All pricing in this brochure is inclusive of GST and are valid until 30 September 2024.

# Green X 12™

THE ADVANCED 4-IN-1 DIGITAL X-RAY IMAGING SYSTEM

# vatech



**5**  
DIFFERENT  
HIGH-RES  
IMAGES

**ALL IN  
1  
SCAN**

**10**  
year  
WARRANTY  
CT SENSOR & TUBE  
\*This promotion is valid for a limited time only.

### SMART FOCUS: 1 SCAN, 5 IMAGES

One 12x8.5, Three 4x4 & One auto pano at the same time!

### 4-IN-1 DIGITAL

Pano | CEPH | 3D CBCT | Model Scan

**MULTI FOV SELECTION**  
**ENDO MODE WITH**  
**0.05mm VOXEL SIZE**

**3d MODEL SCAN**  
**41 LAYERS PANORAMIC**  
**IMAGES: INSIGHT PAN**

**ONE 12X8.5, THREE 4X4**  
**& ONE AUTO PANO**

## EzRay Air Wall

Easy handling and stable positioning



**BEST  
SELLER**

**vatech**  
**EzRay Air™ Wall**

## EzSensor Classic

Your partner for digital



**BEST  
SELLER**

**vatech**  
**EzSensor Classic**

# Green X 16™



**10**  
year  
WARRANTY  
CT SENSOR & TUBE  
\*This promotion is valid for a limited time only.

**Speed, quality, predictability.**  
**NO COMPROMISE**

### The advanced 4-in-1 Digital imaging system

- **49.5µm CMOS Sensor**
- 4 in 1 Digital (Pano/CEPH/CBCT/Model Scan)
- Speed, quality, predictability

### NO COMPROMISE

- Airway volume analysis
- Endodontics module included
- Includes 10 year warranty\*
- Includes PC and monitor

## Own the **vatech** A9™ now...

**UPGRADE** to the Vatech Green X16 or Green X18 within 3 years and we will **BUY BACK** your Vatech A9!\*



## **vatech** A9™

- High quality panoramic
- Integral 8x8 FOV for 3D images
- AI technology with magic pan
- Includes PC and monitor



**10**  
year  
WARRANTY  
CT SENSOR & TUBE  
\*This promotion is valid for a limited time only.

# Green X 16™

- **49.5µm CMOS Sensor**
  - 4 in 1 Digital (Pano/CEPH/CBCT/Model Scan)
  - Speed, quality, predictability
- NO COMPROMISE**
- Airway volume analysis
  - Endodontics module included
  - Includes 10 year warranty\*
  - Includes PC and monitor

\*Terms and conditions apply. Contact us for details.

# 44W of unprecedented power.

## NEW Ti-Max Z Series Air Turbine Handpieces

Z990L/Z890L

Experience like never before, the unprecedented power of the industry's best-in-class 44W air turbine handpiece as it delivers powerful cutting strength.

It makes zirconia removal easy, reducing chairside time and minimising patient discomfort.



## NEW Ti-Max Z2 Series Handpieces

### Z990L

Optic	Model	Order Code	RRP
✓	Z990L	PA23730001	\$1,981

#### Features:

- Head: Standard head
- Power: 44 W\*
- Rotation: 280,000-360,000 min-1
- Head size: ø12.1 x H 12.7 mm
- Water spray: Quattro spray
- Body material: Titanium
- Body coating: DURAGRIP

\*Z990L, Air pressure at 0.42 MPa.



**BUY 2**  
Ti-Max  
Z990L  
and get  
**1 FREE**

### Z890L

Optic	Model	Order Code	RRP
✓	Z890L	PA23720001	\$1,981

#### Features:

- Head: Miniature head
- Power: 41 W
- Rotation: 320,000-400,000 min-1
- Head size: ø10.6 x H 10.8 mm
- Water spray: Quattro spray
- Body material: Titanium
- Body coating: DURAGRIP



**BUY 2**  
Ti-Max  
Z890L  
and get  
**1 FREE**

## NSK Various Combi Pro Bundle

### TWO SYSTEMS IN ONE

Combining ultrasonic scaling and powder therapy for more effective and comfortable oral hygiene and biofilm control.

### Special Offer

#### Includes:

Various Combi Pro Unit

#### PROPHY

- Chamber (supragingival air polishing)
- x 3 handpieces (supra)
- x 4 bottles of FLASH pearl supragingival powder

- PERIO Chamber for subgingival air polishing
- X 2 perio handpieces
- X 2 boxes (40pcs) Sterile disposable nozzle tips
- X 2 160g bottles of Perio-Mate powder

#### SCALING

- X 2 VA2LUX - Optic scaler handpieces
- X 3 scaling tips (G4, G6, G8) and steri case

ICART V to keep everything together!



**VARIOUS COMBI PRO**  
-RRP: \$15,278-  
**\$12,751**  
INCL. GST

**SAVE \$2,536**

# Xceed Limits

## Lara XL

lara XL



**BUY A LARA XL**  
+ GET FREE  
INSTALLATION

**IN STOCK!**

### W&H Lara XL

Eco Dry + technology

High capacity 28L

Unique connectivity and upgradeability

Advanced traceability options

Easy to navigate

Environmentally friendly packaging

### W&H Lisa



Incredible inside and out

Built for highest demands

**EliSense:** Incredibly user friendly

**EliTrace:** Trace and document down to the individual instrument

**EcoDry+:** Eco Dry technology that adapts the drying cycle to the mass of the load. Reducing cycle time, increases the life span of your instruments and optimises energy consumption.



**IN STOCK!**

**BUY A LISA**  
+ TRACEABILITY

**SAVE \$2,000**  
INCL. GST



## W&H Handpiece Maintenance

### W&H Assistina One

assistina ONE

**Fast. Flexible. Thorough.**

Convenient in its handling. Outstanding in its results. With the shortest cycle time in its class and the Quick Connect System, the Assistina One is a top maintenance solution – designed for top dental professionals.

Price:  
**\$2,995** incl GST

**BUY AN ASSISTINA ONE & SAVE \$500**  
INCL. GST



**IN STOCK!**

"The **ONE** that keeps ME safe."